



# Quantitative analysis of residual agricultural chemicals in food by GC-MS/MS

- Quantitative analysis of pesticides in spinach extract -

Product: JMS-TQ4000GC GC-MS/MS System

## Introduction

As "food safety" is recognized as an increasingly important issue on a global scale, many nations have their own regulations on residual agricultural chemicals in food. In Japan, the positive list system, which was enforced at the end of May 2006, stipulates a uniform standard of 10 ppb as a quantity that is considered safe for human health. Under the positive list system, more agricultural chemicals need to be examined, and as a result, techniques capable of accurately and collectively analyzing residual agricultural chemicals in food are in increasing demand. While mass spectrometry (MS) is known for its high detection sensitivity, MS/MS is becoming the mainstream of pesticide analysis for its superior sensitivity and selectivity. The JMS-TQ4000GC, JEOL's latest GC-MS/MS, has a unique ion storage/ejection mechanism within the MS/MS collision cell and incorporates new firmware to support MS/MS analysis with up to 36,000 transitions. In this work, we performed quantitative analysis of residual agricultural chemicals in spinach extract using a JMS-TQ4000GC.

## Experiment

A pesticide standard solution from FUJIFILM Wako Pure Chemical Corporation (PL series) was used that consisted of equal amounts of PL 1, 2, 3, 4, 5, 6, 11, and 12. Afterwards, the solution was diluted to 1, 5, 10, 50 and 100 ppb. PEG 300 was used to protect the pesticides from thermal decomposition in the GC injection liner.

For the sample, 15 g of spinach was processed by using AOAC 2007.01 extraction method, and the resulting extraction solution was mixed with 100 ppb of the standard solution at 9:1. The sample was quantitatively analyzed for 150 pesticides. Table 1 shows the measurement conditions used for the analysis.



GC-MS/MS, JMS-TQ4000GC

Table 1. Measurement conditions								
System	JMS-TQ4000GC (JEOL)							
Ionization mode	EI+: 70eV, 50μΑ							
GC column	VF-5ms, 30m x 0.25mm, 0.25µm							
Oven temp.	50°C (1min)→25°C/min→125°C →10°C/min→300°C							
Inlet temp.	250°C							
Inlet mode	Splitless, 2µL							
He flow	1.0mL/min (Constant Flow)							
MS/MS mode	Peak Dependent SRM							

#### Table 1. Measurement conditions





### Results and discussions

Figure 1 top row shows the data acquired from the original spinach extract while the bottom row shows the data the spinach extraction solution with the pesticides added. MS/MS, with its high mass selectivity, detected agricult chemicals without being affected by contaminants in the spinach extract.

Table 2 shows the quantitative results of 150 pesticides (n=5) and their reproducibility (CV). For most agricultural chemicals, the recovery rate was 70 to 120% and the CV was 10% or less, demonstrating the effectiveness of the . TQ4000GC for pesticide analysis.



Upper: Spinach extract solution only, lower: 10ppb pesticides in the spinach extract solution

Nº Carlor hu			and a second	-	I An
	100	No.	0-0	70	



## Table 2 Quantitative result 150 pesticides in the spinach extract solution (n=5)

No	Compound	PT (min)	Quantitative values (ppb)		)	Ave.			Compound	PT (min)	Quantitative values (ppb)					Ave.	01/1%)		
100.	Compound	KT (mm)	No.1	No.2	No.3	No.4	No.5	(ppb)	CV (%)	NO.	Compound	KI (min)	No.1	No.2	No.3	No.4	No.5	(ppb)	CV (%)
1	Atrazine	12.37	10.0	10.0	9.7	10.4	9.2	9.9	4.4	76	Pyriminobac methyl 2	12.88	10.8	10.9	11.0	10.4	10.8	10.8	2.1
2	Benthiocarb	14.68	10.1	10.7	11.6	11.5	10.8	10.9	5.8	77	Pyroquilon	19.53	10.3	9.2	9.7	9.7	10.0	9.8	4.2
3	cis-Permethrin	20.68	12.5	12.8	12.8	12.3	12.6	12.6	1.9	78	Tetradifon	13.87	7.5	9.3	9.4	8.5	8.1	8.6	9.1
4	Diflufenican	18.31	11.7	11.6	12.0	11.2	10.7	11.4	4.3	79	Tolclofos-methyl	16.61	9.2	9.2	9.2	9.2	9.1	9.2	0.6
5	Fenamiphos	16.24	9.1	11.6	10.4	10.9	9.5	10.3	9.8	80	Uniconazole P	13.66	11.5	12.4	10.8	11.5	10.6	11.4	6.3
6	Fenarimol	20.20	11.1	11.1	11.1	11.3	11.1	11.1	0.7	81	Acetochlor	17.83	9.3	8.7	8.8	9.3	8.9	9.0	3.3
7	Fenpropimorph	14.74	11.0	11.4	10.7	11.8	10.6	11.1	4.6	82	Benalaxyl	13.57	10.7	9.9	10.3	10.2	9.8	10.2	3.5
8	Norflurazon	17.94	10.7	10.2	10.6	10.2	10.7	10.5	2.3	83	Benfuresate	11.71	9.6	9.3	9.4	9.7	9.6	9.5	1.6
9	Oxadiazon	16.50	10.6	9.9	11.0	10.8	10.3	10.5	4.0	84	Cadusafos	11.37	8.2	8.3	8.0	8.6	7.9	8.2	3.5
10	Penconazole	15.45	11.7	10.8	10.9	11.3	11.2	11.2	3.0	85	Chlorpropham	15.61	8.9	8.5	8.2	9.1	8.9	8.7	4.4
11	Pendimethalin	15.29	12.9	11.7	12.3	11.9	12.5	12.3	3.7	86	Diclocymet 1	15.91	10.4	10.8	10.4	10.9	10.2	10.5	2.9
12	Procymidone	15.70	10.6	10.8	12.4	9.4	9.7	10.6	11.2	87	Diclocymet 2	15.40	10.0	10.4	10.8	9.9	10.0	10.2	3.4
13	Spiroxamine1	13.78	8.6	8.3	8.2	8.6	8.5	8.4	2.1	88	Dimethametryn	14.49	12.0	10.7	11.7	11.6	11.4	11.5	4.3
14	Spiroxamine2	14.31	9.2	8.7	9.0	9.4	8.5	9.0	3.9	89	Esprocarb	21.83	9.3	9.4	9.3	9.7	9.1	9.4	2.5
15	Tefluthrine	12.92	10.2	9.9	10.3	10.5	10.1	10.2	2.2	90	Etofenprox	16.08	12.8	12.6	12.7	12.9	12.4	12.7	1.5
16	Terbufos	12.63	8.3	7.9	7.9	8.7	7.6	8.1	5.0	91	Fenothiocarb	13.27	10.4	10.3	10.4	10.7	11.0	10.6	2.9
17	Terbutryn	14.31	10.5	9.6	10.1	10.4	9.7	10.0	4.0	92	Iprobenfos	16.92	9.6	9.0	10.0	10.0	9.2	9.6	4.7
18	trans-Permethrin	20.68	12.5	12.8	12.8	12.3	12.6	12.6	1.9	93	Isoxathion	17.61	6.6	8.0	8.0	6.7	6.9	7.3	9.8
19	Alachlor	13.84	9.7	9.0	8.9	9.3	8.7	9.1	4.5	94	Mepronil	14.04	11.4	11.1	11.8	11.6	11.5	11.5	2.2
20	Buprofezin	16.68	11.2	11.7	11.1	11.9	11.9	11.6	3.2	95	Prometryn	10.85	10.2	9.7	10.0	10.9	9.8	10.1	4.8
21	cis-chiortenvinphos	15.49	0.3	/.6	/.9	/.6	0.4	/.2	10.4	96	Propachior	16.41	9.4	/.6	/.0	8.0	/.5	7.9	11.5
22	cyproconazole 1	17.03	12.1	12.4	11.2	9.9	12.2	11.6	9.1	97	Prothiotos	15.99	8.7	9.7	8.9	9.3	9.2	9.1	4.4
23	Cyproconazole 2	17.03	11.6	11.0	12.0	12.3	11.3	11.7	4.6	98	Pyritenox1	12.89	10.2	10.0	10.6	10.5	10.6	10.4	2.3
24	pirenoconazole 1	23.10	11.3	11.2	11.3	11.0	11.3	11.2	1.4	99	Pyrimetnanil	13.94	9.4	9.4	9.6	10.5	8.9	9.6	6.0
25	Ditenoconazole 2	23.18	11.7	11.1	11.2	11.3	11.4	11.3	2.0	100	simetryn	13.06	11.0	10.6	9.4	10.2	10.7	10.4	6.0
26	Ethion	17.34	10.2	10.6	10.9	10.0	9.9	10.3	4.0	101	I erbacil	14.82	8.7	10.0	9.9	10.2	9.2	9.6	6.6
27	renitrothion	14.33	8.3	9.0	10.7	9.3	9.1	9.3	9.5	102	l etraconazole	18.31	11.9	11.6	11.7	10.9	13.4	11.9	7.7
28	Fenthion	14.73	8.9	10.0	8.7	9.4	8.4	9.1	7.1	103	I nenyionior	16.58	/.1	8.1	8.8	8.0	/.8	8.0	/.3
29	Fluridone	22.18	11.1	12.0	12.0	11.7	11.6	11.7	3.3	104		16.66	9.0	10.4	9.7	9.8	8.3	9.4	8.4
30	Hexazinone	18.24	9.5	9.9	10.0	9.5	9.0	9.6	3.9	105	Tricyclazole	15.72	4.4	6.9	6.3	5.4	4.7	5.5	18.8
31	Isotenphos oxon	14.79	10.1	12.8	13.3	10.4	10.6	11.5	12.9	106	Zoxamide (decomposed)	13.40	10.5	11.1	11.0	10.3	9.6	10.5	5.7
32	Isophenphos	15.43	9.7	8.5	10.0	9.8	10.1	9.6	6.5	107	Benoxacor	14.44	4.2	4.7	5.0	4.5	3.8	4.4	10.1
33	Isoprothiolane	16.44	11.1	11.2	11.0	10.1	12.4	11.2	7.4	108	Bromacil	13.70	6.6	8.1	8.1	8.4	7.9	7.8	9.0
34	Propargite 1	18.31	6.7	7.2	7.4	7.0	7.3	7.1	4.1	109	Bromobutide	16.17	11.7	9.7	10.5	11.8	10.2	10.8	8.7
35	Propargite 2	18.31	6.9	7.3	7.4	7.2	6.8	7.1	3.9	110	Butamifos	13.56	10.3	12.0	11.7	11.0	10.9	11.2	6.3
36	Propiconazole 1	17.93	12.6	12.9	12.2	12.1	12.0	12.4	3.1	111	Dichlofenthion	15.12	9.7	9.3	9.5	10.2	9.1	9.6	4.6
37	Propiconazole 2	18.07	13.5	12.4	12.3	12.9	13.0	12.8	4.0	112	Diphenamid	16.42	12.2	11.5	12.2	11.9	12.1	12.0	2.3
38	Propyzamide	12.73	11.0	10.6	11.1	10.4	10.4	10.7	3.2	113	Hexaconazole	14.00	10.0	11.6	11.4	9.9	11.4	10.9	7.7
39	Pyriproxyten	19.72	11.0	10.9	10.9	11.1	10.7	10.9	1.3	114	Merenoxam	16.33	16.3	15.5	17.4	16.6	16.1	16.4	4.4
40	trans-Chiorfenvinphos	15.25	7.6	9.4	8.4	8.5	6.4	8.1	14.1	115	Napropamide	17.40	10.6	10.6	10.8	10.5	10.1	10.5	2.5
41	Iriadimenol 1	15.70	12.2	10.0	10.5	10.5	10.5	10.7	8.1	116	Oxadixyi	16.07	9.2	10.9	10.2	9.6	9.1	9.8	7.5
42	Iriadimenol 2	15.84	12.8	10.6	11.9	10.7	11.0	11.4	8.1	11/	Paciobutrazoi	19.31	11.8	12.0	12.8	12.3	12.1	12.2	3.1
43		13.13	7.4	7.9	8.1	8.3	7.4	7.8	5.0	118	Phenothrin 1	19.42	13.8	10.0	11.1	12.0	10.5	11.5	13.3
44	Vinciozoline	13.76	7.2	8.0	8.3	7.5	7.6	10.0	5.6	119	Phenothrin 2	18.98	11.5	11.2	11.0	11.0	11.2	11.2	1.8
45	Acetamiprio	18.92	10.3	10.9	10.7	10.3	11.5	10.8	4.8	120	Piperopnos Dia hudia i agrada 1	12.94	9.9	10.0	11.3	10.9	10.5	10.5	5.5
46	Allethrin 1	15.43	10.0	9.5	11.3	10.9	11.7	10.7	8.5	121	Pronydrojasmon 1	13.23	8.8	8.5	8.4	8.9	8.1	8.5	3.8
4/	Allethrin 2	15.51	9.0	10.5	9.2	9.1	10.3	9.6	1.1	122	Pronydrojasmon 2	12.42	11.8	12.2	11.2	14.5	12.2	12.4	10.3
48	Bitertanol 1	20.00	14.1	14.1	14.0	14.1	14.0	14.0	0.5	123	Propazine	15.03	9.7	8.5	9.3	9.8	8.8	9.2	0.0
49	Promopropulato	20.76	10.8	10.5	10.7	10.2	11.2	10.2	3.0	124	Pyributical D	15.51	10.6	12.7	11.1	10.0	12.9	10.2	5.5
50	Chlorobonzilsto	17.00	2.0	10.0	10.7	10.3	7.8 0 C	10.2	4./	125	Cuinalphos	10.14	10.0	7.9	11.1	11.0	10.1	11.0	7.5
51	Chlorovrifos	17.19	8./	9.4	10.1	9.4	8.0	9.3	0.9	126	Quinalphos	19.10	10.1	10.7	11.8	10.0	10.1	11.0	1.3
52	Overfluorfon	16 50	0.7	7.4 12.0	0.0	3.0	0.0	12.0	4.0	12/	Tolfonpurad	24.03	11.1	11.7	11.4	11.5	11.4	11 5	2./
53	Darathion	14 00	12.2	14.4	13.0	12.4	12.0	12.0	5.0	128	Aldrin	16.76	7 /	211.3	9.7	211.3	9.2	27	10.5
55	Piriminhos mothyl	14.80	13.2	14.4 9 0	9.5	13.0	12.0	9.6	3.0	129	cis-Chlordano	1/1 07	7.4 g 5	0.5	9.6	0.5	9.2	0./ Q 1	57
56	Propanil	13.67	10.1	10.1	10.6	10.1	10.5	10.2	1.0	101	Dicofol	16 70	11 7	10.7	10.6	10.0	10 0	11.0	2.0
57	Dyridahan	20.07	10.1	12.0	12.0	11 5	11 1	11.0	2.5	101	Dieldrin	17 10	211.7	8.2	10.0	7 4	10.0	2 0 11.0	1/1 /
50	Quinoxyfen	18.00	11 7	11.5	12.1	11.7	11 5	11.0	2.7	122	Endrin	14.07	0.2	10.0	10.4	10.6	20.0	9.0	10.7
38	Simazino	12.00	11./	10.4	12.1	10.2	10.2	11./	2.1	133	Hentachlor	16.02	5.0	10.0	10.7	10.0	6.0	5./	67
23	Tohuconazolo	10.20	7.3	10.4	7.0	11.1	10.2	3.9	4.9	104	trans-Chlordano	14.40	9.5	9.0	9.7	9.7	9.0	9.5	2.0
61	Triadimofon	14.58	11.2	10.0	10.5	10.5	10.3	10.5	3.0	135	1 Naphthylasstamid-	14.40	2.2	12.0	3.7	12 6	1/1 1	12.0	2.6
62	Triazonhos	17.65	7 4	10.4	10.1	10.0	7.0	10.3	7.0	130	Promonhos other	16.00	14.4	13.2	14.3	12:0	0 4.1	12:2	3.0
62	Amotorn	14.00	10.7	9.1	8.5	8.0	10.4	8.U	9.9	137	Carbovin	15.00	8.9	8.0	9.5	<u> </u>	8.0	8.9	4.1
03	Aracoparala	16.00	12.0	10.1	14.5	10.1	12.2	12.0	2.2	138	Chlorbonsido	15.98	3.5	5.0	5.4	3./	3.3	3.5	1.0
65	Pupirimato	16.65	11.0	10.4	11.2	11 1	10.4	10.0	2.2	139	Chlorofenson	2/ 00	10.0	2.7	7.3 8.7	10.0	7.0	7.0	10.5
65	pupifimate Putachlor	16.00	11.0	10.4	10.5	0.7	10.4	10.8	5.8	140	chorosono	24.99	0.8	0.0	0./	0.2	7.4	7.9	57
67	Chlorthal dimethyl	14.72	5.9	7.ŏ	10.0	5./	9.4	3.3	4.4	141	Disulfoton	12.47	0.0	0.3	0.9	0.9	7.9	0.4	3.7
67	Children dimethyl	12.10	9.4	8.9	9.4	0.8	8.5	9.0	4.4	142	Disulloton	10.01	1.5	1.4	11.0	0.1	11.3	11.4	4.0
08	Diothofoncort	14.02	12.2	8.3	8.4	10.0	10.3	/.5	12.1	143	Epoxiconazole	14.07	10.2	10.0	10.2	10.5	11.3	10.0	2.1
70	Dimoningrate	15.74	12.3	11.8	12.6	10.9	12.4	12.2	0.0	144	Echorumesate	14.3/	11.2	11.1	10.2	10.5	9.3	11.4	4.0
70	Dimepiperate	13.74	10.0	11.4	10.1	12.4	13.4	12.8	7.0	145	Flutriafol	15.17	12.6	11.1	11.7	11.2	11.1	11.4	4./
71	primetrienamiù Etovazolo	10.02	11.0	5.0	11.1	10.4	3.5	3.3	5.0	140	Isazanhas	12.06	10.1	11.3	11./	10.0	11.5	11.0	4./
72	Elucation	17.03	11.0	10.0	11.1	11.2	10.0	11.0	0.9	14/	ISazophos Rhorato	11.70	10.1	9.4	9.0	10.4	9.4	9.8	4.0
73	Fiuaci ypyrim	10.14	10.0	10.8	11.0	11.2	10.0	11.0	2.3	148	Priorate	10.07	0.5	8.5	11.0	0.5	0.4	8.Z	4.9
74	LendCli Duriminghager attacks	17.00	10.0	12.3	11.4	12.7	10.0	12.5	0.0	149	Propaphos	15.00	10.9	11.5	11.8	10.7	10.0	11.1	0.1
/5	Pyriminobac metnyi 1	11.90	12.8	12.3	12.5	12.7	12.0	12.0	1.5	120	Fropaprios	13.83	10.8	11.2	11./	10.7	10.9	11.1	5.5

11 Dearborn Road, Peabody, MA 01960 Tel: (978) 535-5900 • Fax: (978) 536-2205 ms@jeol.com • www.jeolusa.com