


Application of ^{14}N decoupled HETCOR for ^{19}F -mediated amide hydrogen bonds with ROYALPROBE™ HFX

Product used : Nuclear Magnetic Resonance (NMR)

Fluorine and nitrogen containing organic compounds are often found in a variety of fields such as pharmaceutical and material sciences. Such compounds show unique properties based on their conformations and orientations, which are often induced by hydrogen bonds between N-H and F atoms. Therefore, analysis of molecular conformations and inter-molecular interactions are important. Here, we report ^{14}N decoupled ^1H - ^{19}F HETCOR as an example of hydrogen bonds analysis between N-H and F for fluorinated benzanilide (**1**) as a demonstration.



fluorinated benzanilide (**1**)

Measurements of ^1H NMR with ^{14}N decoupling

NMR spectroscopy is a powerful tool for analysis of hydrogen bonds, and its J coupling value provides the information of hydrogen bond. Fig.1(A) shows ^1H NMR spectrum of **1**. Amide ^1H signal is observed in the lowest magnetic field region, ca 8.77 ppm. In this case, severe line broadening is arising from quadrupolar ^{14}N . In contrast, measurement of ^1H NMR under the ^{14}N decoupling condition, resulted in the sharp amide ^1H signal observation as double of a doublet (Fig.1(B)). This splitting signal is derived from bifurcated hydrogen bond bridges between N-H (11) and two F (F5, F6) atoms.




Fig.1 (A) ^1H NMR spectrum of **1** (B) $^1\text{H}\{^{14}\text{N}\}$ NMR spectrum of **1**

Sample: 36 mg fluorinated benzanilide in chloroform- d

Equipment: JNM-ECZ400S with ROYALPROBE™ HFX*

*Special modification for ^{14}N nuclei required.

Measurements of ^1H - ^{19}F HETCOR with ^{14}N decoupling

In order to confirm two hydrogen bonds between the amide ^1H and ^{19}F atoms, we collected ^1H - ^{19}F HETCOR under the ^{14}N decoupling condition. 1D-HETCOR spectra of **1** are shown in Fig.2. As expected, we could observe cross peaks to amide ^1H with both F atoms (F5, F6). Furthermore, by using anti phase signals in 1D-HETCOR spectra, J constants were clearly observed as ca. 16 Hz (amide-F5) and 3 Hz (amide-F6), respectively in the case of ^{14}N decoupling conditions. 2D-HETCOR spectrum with ^{14}N decoupling is also available, which is shown in Fig.3.




Fig.2 1D-HETCOR spectra of **1** with or without ^{14}N decoupling
(A: selective excitation at F5; B: selective excitation at F6).

Scans: 4




Fig.3 ^{14}N decoupling 2D-HETCOR spectrum of **1**

Reference

G. N. Manjunatha Reddy, M. V. Vasantha Kumar, T. N. Guru Row, N. Suryaprakash, *Phys. Chem. Chem. Phys.*, **12**, 13232–13237 (2010).

Copyright © 2019 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

